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**A**bridge dividing attachment made out of ZrO<sub>2</sub> paves the way for otherwise hard to solve solutions. Whether with the mass-produced cercon-link attachment or, as it is in this case, an individually cut pivot - there are a lot of very nice and simple solutions to be found and designed. It has been shown that an allform extension at the basal area of the pivot affects the static of the adhesive joint positively. It must be pointed out that scientific tests and examinations are only available in connection with the cercon-link attachment.

empfohlene Bohrer: Zylinder 882-014TC, 881-014TC, 882-014F, 881-014F, 899-027C, 899-027F aus ökoDENT Präpset 2

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You should have a look at picture 3 under the premise that currently maximal two bridge elements for spanning a gap are allowed. Do not tell me that dental technicians are not daring. This all-ceramic bridge produced with three parts, which were connected with a ceramic adhesive, fired and faced later on, staid in the patient's mouth for 18 months. At this point a further note to picture 3: These delicate margins demonstrate part of the opportunities provided to you by the material ZrO<sub>2</sub>. The enormous homogeneity of this material paves the all-ceramic way to unprecedented finest crown margins.


empfohlene Bohrer: Zylinder 882-014TC, 881-014TC, 882-014F, 881-014F, 899-027C, 899-027F aus ökoDENT Präziset 2

## Bildbeschreibung / captures

- 1 ZrO<sub>2</sub>Abutment / ZrO<sub>2</sub>-abutment
- 2 Cercon-Gerüst vorbereitet zum Verblenden  
Cercon frame prepared for facing
- 3 fertiggestellte Arbeit basal / finished work from basal
- 4 fertiggestellte Arbeit okklusal / finished work from occlusal
- 5 okklusale Ansicht in situ / occlusal view in situ

# inlay-bridges discover the possibilities

1

2

3

4

5

Autor: Dentist Dr. Jörg Jürgens

**D**iscover the possibilities – an often-cited sentence which hits the nail on the head! The in wax prepared inlay 37 with individual intracoronal matrix subsequently bears the integrated male part of attachment, which is a restoration of the hermetic tooth 36 also made out of ZrO<sub>2</sub>.

In this case the unstable base was preferred to a fixed blocking. From basal in the interdentially approximal area one can identify the extremely homogeneous surface which can be achieved by diamonds of the finest grain with the help of a water-cooled turbine. This is also an advantage with little available space in the area of the connector.

It is to be mentioned that depending on the individual wish of the patient, dentist, or dental technician the form of the marginal areas and, if necessary, also the complete facing can be achieved under utilisation of the Cercon-Ceram-Press-Ceramic attuned to the cercon system. Aesthetically a complete success!

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

## Do you want a little bit more?

Autor: Dentist Dr. Angela Langschwager

**D**o you want to get more for your money? All-ceramic locking bar constructions, as it is in this case, are still likely to be an exceptional technology at that point of time. However, a prospect of the technical realisation of such works is to be exemplified briefly:

16/17 receives a provision on the basis of a telescope

13 all-ceramic crown with distal locking bar pivot

12/22 all-ceramic bridge elements

11.21.23 all-ceramic crowns

Before "it comes to the crunch concerning the locking bar" a preliminary ceramic complete insertion is incessant because the bridge must not be treated thermally after the electroplating process of the distal pivot at 13. The now following production of the zirconium oxide case and blade prove to be a very uncomplicated process as it can be read in "Tipps und Tricks". Less is yet more as it is generally known! In order to reduce the abrasion in the area of the locking bar pivot the galvanised part of the primary frame was galvanised once more in this connection. This fitting part was adhered to the already faced secondary frame. Consequently you will receive a comfortable attachment fitting.

empfohlene Bohrer: Zylinder 882-0147C, 881-0147C, 882-014F, 881-014F, 899-027C, 899-027F aus ökoDENT Präpset 2

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Zirconium frames excuse discoloured stumps or pin constructions made of noble metal even under a wall thickness of 0.3 to 0.4mm. Even with by all means present translucency of the frame material it blocks very reliable present discoloured irritations with the help of the base. As it is present in this case, an implant abutment can be integrated optimally. As in nearly all cases the cementation is carried out the conventional way. In contrast to conventional all-ceramics the one at hand allows due to the enormous bending strength a trial wearing and the grinding of possible premature contacts before the definite fixation. A control of the margins is possible with the help of the radiographic opacity of the zirconium oxide.

empfohlene Bohrer: Zylinder 882-014C, 881-014C, 882-014F, 881-014F, 899-027C, 899-027F aus ökoDENT Präpset 2

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# Cercon/Solaris clever bolt-technology

Autor: Dentist Dr. A Langschwager

Two technical, highly interesting highlights found their way into this work with the remaining teeth 13, 33, 23, 17 and 27. On the one hand the problem had to be solved to reduce the abrasion of a noble metal circulation at the anterior primary frame made out of zirconium oxide and on the other hand to find a way to activate the locking bar blade comfortably if necessary.

A firing on of the ceramic Cercon-Ceram 5 in the area of the circulation and the subsequent re-cutting proved to be successful in this connection. An abrasion could apparently not be found any more. The realisation of a locking bar blade, which can be activated, was also easy to convert. The technician opened the cut frame still being chalky in the area of the locking bar pivot from occlusive with the help of a spiral bone burr.

After the sintering procedure a finishing of the circulations and the pivots followed with the help of a water-cooled turbine in the cutting device.

Afterwards also under water-cooling the technician opened up manually the intakes for the lobes of the locking bar blades with the help of a special tool equipped with diamonds.

After a ceramic complete trial the electroplating of the inserted friction-retained pins was carried out.

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
**A**ll-ceramic version or conventional provision with a telescope? There are two options to gain a solution in the present case.

1. The conventional version with a circular resin-faced bridge on five telescoped primary crown out of non-noble metal.
  2. The all-ceramic solution with a galvanised cap which was adhered to the faced zirconium frame.
- The solutions will be introduced on the following double page.

# telescop plus cercon?

## 1 case – 2 solutions

Autor: Dentist Dr. Antonia Fuchs

There are two options for a solution in the case at hand. For the conventional version the technician made a circular resin-faced bridge on five telescoped primary non-noble metal crowns.

For the all-ceramic solution a secondary galvanised cap was used which was adhered afterwards to the faced zirconium frame.

If you focus only on the point "friction", the decision for one or the other version is not hard to make. Even if the decision was made for a telescoped solution as it was the case with this patient, the design of the primary crowns with circular 20-cutting eased it up for you. In this connection three points are of significant importance for the long-run success of such restorations:

1. The creation of a gap as small as possible between the primary and secondary crowns (<5mm) can be achieved with a singular spraying out of the conducting silver varnish with the help of an airbrush.

2. The definite insertion (cementation) of the primary crowns and the adherence of the galvanised secondary crowns in the tertiary frame in the mouth at the same time. This will eliminate all possible sources of error of the master model and the impression (in connection with the not excludable minor distortions of the used materials).

3. Important! Aim for a slight "clearance fit" between the tertiary frame and the secondary crowns which allows for compensating the finest dimension changes of the impressions and the models in the mouth. Do you have already "experienced" the action of a galvanised cap on an all-ceramic primary crown?! With this procedure you can achieve to design a defined friction/adhesion reproducible. You spare yourself a time-consuming fitting on and adjustment of the secondary fit and you will deploy all advantages of this innovative material. Even the topic of cold welding as you experience it occasionally in the conventional area between primary and secondary crowns will be "no longer an issue" with the usage of all-ceramic primary crowns.

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**T**hese galvanised secondary parts were produced by the laboratory on four all-ceramic 2° cut zirconium primary parts. The surfaces exposed by the primary crowns when processed in a cutting device with water-cooled turbine and decent diamond set will "whet your appetite". You will achieve a "nacreous" mirror finish. A circular shoulder with rounded inner edge on the one hand will increase the stability of the secondary crowns and on the other hand will create the base for a comfortable support of the tertiary construction. Due to the extension of the blank sizes around the Cercon-base 47 (47mm anatomic length) it would be possible to produce the tertiary frame as one piece.

In this case the bridge parts from 33 - 35 and from 43 - 45 were cut and sintered with one individual pivot each mesio-proximal of the canine tooth crowns. The subsequently modelled and cut segment 32 - 42 was adhered to the parallel pivot with a ceramic adhesive by the technician. The compatible adherence of the three zirconium segments allowed for a cosmetically invisible covering of the basal areas with the facing ceramic CerconCeram5 in the course of the facing.

empfohlene Bohrer: Zylinder B82-014TC, B81-014TC, B82-014F, B81-014F, B99-027C, B99-027F aus akadENT Präpset 2

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Note: Such a ceramic adhesive was not yet tested in composition of a zirconium restoration and no otherwise test results in this connection are available. Furthermore, it has to be mentioned that nowadays only two intermediate elements for spanning a gap are allowed.

It is quite encouraging that a good many of such works are required by patients and dentists and that those are functioning by all means. The inventive raring to go in the laboratories and the noticeable momentum resulting from it are motivation and inspiration to me.

Time and again during our technical instructions in the laboratories we experience that the technicians are not interested in working within the areas of the limiting values with these high-tech products. In fact it seems that stimulation exists for knowing the limiting values and going beyond them. An example: I often experience that a training model of a 3-element bridge causes weary smiles at the most from the "all-ceramic enthusiastic" technicians and nobody is really "keeled over" by it. Later on the models which are presented to me with a winning smile by the other side are tough ones, such as e. g. Maryland Bridges and Inlay Bridges, locking bar works, implant borne structures, attachment technique etc.

And now we really have to ask ourselves where we would stand medical-technical wise if nobody has dared to go beyond the limits time and again. What happens on the dentistry market at the moment inspires and encourages.

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

For the removal of conventional crowns and bridges a so-called crown cutter is used in general. Here it concerns tungsten-carbide burs which rip up metal alloys, titanium, as well as facings out of low-fusing ceramic and polymer preferably with 60,000 to 80,000 rpm (revolutions per minute); whereas they are unsuitable for all-ceramic out of extremely hard ceramic, e.g. zirconium oxide.

By today's view ZrO<sub>2</sub> works are ripped up ideally with the help of a coarse-grained diamond grinder (150 green) in a high-speed grinder with 40,000 rpm under ample water-cooling. This ripping up process seems to cause a higher tool wear and appeared at first a bit more difficult to the untrained user but from our point of view it avoids a pulpal damage. The fact is worth mentioning that after ripping up the ZrO<sub>2</sub> construction the bending up procedure for the separation into parts can be carried out often without noteworthy ceramic splintering (see crown 11). A diamond grinder too coarse (ISO black) and/or a cutting procedure at a speed too high does not only increase the risk for an irreversible pulpal damage but also prevents the above mentioned phenomenon (see crown 21).

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