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



Basically, every implant-prosthetic rehabilitation is a special challenge for the prosthetist. On the one hand side, heavy angular inclinations and divergencies have to be compensated for, on the other hand side, the ideal material must be selected. Lastly, superstructures either have to be processed in the minimum layer strength or must bridge enormous distances.

Zirconium oxide ceramics is very suitable for the whole implant prosthetic due to its material properties, its biocompatibility and its outstanding aesthetics (see page 24, prosthetic treatment Dr. Arvid Langschwager and ZM Präzisionsdentaltechnik GmbH). Heavy angular inclinations can be levelled by individually milled abutments from zirconium oxide (see page 24).

These ceramic abutments are by far superior to the individual abutments produced with alloys. The surface can be polished up to 3µm and the abutments are not subject to static charges as metals are. Therefore, adhesion of dental plaque and bacteria can be avoided as far as possible. (see page 24). Furthermore, the enormous stability of the material provides for an extremely thin material layer, at the same time providing for almost any angle possible. The white colour of the abutment leads to a clearly defined lucency of the gingival garland (see page 24).

With this method, complicated cases can be solved in an aesthetic way as well. Furthermore, the enormous tensile strength of the material provides for fully tensioned telescopic bridge constructions from complete ceramics. These pieces are durable, are not subject to modifications due to attrition or discolouration and have a low allergic potential due to the few materials used (see "FACE and more", special edition 1/2005). Moreover, the often quoted argument that it is impossible at the moment to mill complex structures or great distances from zirconium oxide, is totally unfounded and is based largely on the lack of appropriate knowledge.

Of course, the size of the blank pieces limits the dimensions of the structure, but a separating bolder and a visionary adhesion technique before the ceramic baking make it possible to combine any elements, so that a limitation for pieces from zirconium oxide virtually does not exist anymore (see page 8 and page 36 of the main edition I "FACE and more"). Therefore, it is quite possible to realise circular millings, bolts (see page 44/45 of the main edition I "FACE and more") and enormous dimensions. Furthermore, it is possible to combine frameworks of different brightnesses and colourings (yes, zirconium oxide can be dyed!) with each other in order to avoid opaque baking (see page 6/7 prosthetic treatment Dr. Arvid Langschwager and ZM Präzisionsdentaltechnik GmbH).

Particular importance is attached to the fact that the adhesive is a temperature-resistant adhesive, which "burns" the pieces into each other as it were. The advantages of this technology are obvious. The costs, which increase enormously with the size of the blank pieces, are minimised. It is much more cost effective to use several smaller pieces to full capacity than to use oversized blank pieces. In addition, the individual parts can be better placed in the blank piece, which results in a tremendous saving of blank pieces.

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


fullceramic total supply

In the case covered here, it would not have been possible to mill the framework with any system available on the market to date (see page 8 prosthetic treatment Dr. Arvid Langschwager und ZM Präzisionsdentaltechnik GmbH).

Due to this limitation, the wrong compromises are made very often and, for instance, the posts of the branding ceramics are not sufficiently supported by the framework – this leads to failures, which are wrongly attributed to an insufficient connection with the framework. However, a particular material can only be used as a standard in dentistry, when it is not limited to only a few indications, but when it can be used for almost any form of dental treatment!

As innovative and progressive prosthetist and dentist, one is obliged to offer the best (state of the art) treatment according to the state of the science to the patient, more so if clear advantages can be achieved for the patient with regard to material-technical (stability/aesthetics) and health aspects (plaque accumulation allergies) by the choice of material.

However, the aim is to offer these advantages to the patient on a timely basis without having to wait until the sometimes very slow industry is ready to adjust the systems to the requirements of everyday life. In this connection, the question is raised how the profitability will develop with more complex machines and attention has to be paid so that the added value remains in the hands of the dental technicians and dentists, so that they are not subject to a doctrine concerning price and choice of material sometime in the future.

Very good reasons, why this visionary adhesive technology before the ceramic baking is what all users of zirconium oxide have been waiting for for a long time. Of course, many tests and some attempts have been made to develop a temperature-resistant adhesive for zirconium oxide. However, the baking temperature for screening high-quality ceramics could not be reached to start with, time and again bubbles and cracks formed, let alone fractures.

In cooperation with DCM GmbH, Dr. Arvid Langschwager and the working group Zirconium Oxide of ZM Präzisionsdentaltechnik GmbH (see main edition I "FACE and more"), we can pride ourselves with a long-term experience of more than three years. The partially very complex structures with both interdental and occlusal adhesion (see page 18/19 prosthetic treatment Dr. Arvid Langschwager and ZM Präzisionsdentaltechnik GmbH) were tested and all cuttings and grinding samples bore up against a microscopic examination. We do apply this technology to our patients as a routine treatment.

For one year, the following four dental laboratories have been working successfully with this adhesive technology:

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We (practice-clinic Kielhorn + Kielhorn) are pioneers in the use of zirconium oxide. However, after having realised the advantages of the material, we learned the limits very rapidly as well, despite of permanent innovations (we already introduced the implant-borne, milled zirconium oxide bridge to the technical press three years ago). Moreover, we discovered unfortunately that the industry tinkers around on its own without involving pragmatists and professional people.

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Dental Design Reichert GmbH