ZrO₂ – geklebte Abutments

ZrO₂ – bonded abutments / ZrO₂ – abutments

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XIVE ZrO₂ - Abutments that have been coated using **CERAM PRESS** allow for a good alignment with regard to the course of the gingiva and the functionality of the conical abutment. Due to its quick production mode, this method is very famous. Unfortunately, the porosities of the over-pressed ceramics and the material's brittleness do not allow for a perfect surface refinement that should serve as prerequisite for an over-galvanisation. This material characteristic considerably complicates further processing.

Abutments that are bonded with a conical coping allow for an individual design, as well. By the „hot bonding“ a homogeneous and stable bonding has been created between coping and structure. This bonding allows for an unproblematic handling concerning the design and the pre-conditioning for the over-galvanisation. Even a coping that is ground through does still guarantee good stability and good superfinishing of the ZrO₂ surface (see fig. 1,2,3 below).

This is a little bit more time-consuming, however, absolutely more worthwhile.

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www.dcm-management.de

Profis schulen Profis
"Professionals coach professionals"

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zircon hotglue - technology

The DCM Competence Team "Zirconium Dioxide" has been successful in developing a paste for the hot pasting of zirconium dioxide before the ceramic fire (zircon hot paste).

In the tests, in particular with laboratories outside of the competence team, it has become clear that only the use of the paste according to a steady technology ensures the successful production of multi-limbed inter-dental and occlusally pasted bridges, as well as ZrO₂ bars and abutments.

Therefore it is necessary from our point of view in order to secure the quality to train the laboratories which want to use the paste in this technology.

This training includes the following working steps for you

Pic. 1: Finished model
ZrO₂ 11-13 bridge for the occlusal extension primary
14 inter-dental ceramic module
23 Xive ZrO₂-Abutment for overpressing with
ceram press
26 Xive ZrO₂-Abutment for 2° conical crown
primarily for heat pasting

Pic. 4: secondary parts to be made during the training
sessions (1st day)

Pic. 5: Adjustment of the overpressing and the secondary
parts (2nd day)

Pic. 6: Completion of the end product after heat pasting
and post-processing through water cooled turbines






The fees include the following training material:

- 1 Model
- 1 Bridge (cerconbase)
- 2 Xive ZrO₂-Abutments
- Cercon ceram press pellets
- Cercon base blanks

Turbines, milling devices and paste are made available for the duration of the training session.

In order to retain a high degree of quality during the training session, the training sessions will not exceed a maximum number of 10 participants.

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Over-pressed abutments vs. bonded ZrO₂ abutments

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Over-pressed, compact abutments with a slightly dull surface, thus very susceptible to plaque

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tiny, bonded abutments with a homogeneous, polished surface, thus plaque-resisting

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